

TANNAS CO.

4800 JAMES SAVAGE ROAD
MIDLAND, MI 48642
USA



"Tomorrow's Instruments... Today"™

SPECIFICATIONS

Scanning Brookfield PlusFour (SB+4) Bath

Principle:

- Liquid bath for analyzing low temperature viscosity and Pumpability/Gelation Index on fluids and lubricants.
- Cost-effective approach for labs with minimal bench space and low to moderate testing needs.

Dimensions:

- Bench-top: ~125 lbs. (~28"(w) x 20"(d) x 34"(h))

Cooling Capabilities:

- Approximately +30°C to -45°C
- Programmable cooling rate: 0.1° to 10°C/hour

Bath Control:

- Programmable: $\pm 0.1^{\circ}\text{C}$

Bath Medium:

- Methanol is recommended.

Safety:

- High temperature Cut-out
- Low Liquid Level Cut-out
- Multiple fuses

Standard Voltage:

- 120 VAC, 60 Hz., single-phase *(Also available in 220 VAC and 50 Hz.)*

Test Methods/Specifications:

- ASTM D 5133, new D 7110 method for used and soot-containing oils
- Required for ILSAC GF-2, GF-3, GF-4; API 'SL', 'SM' category
- Recommended for SAE J300 Viscosity Classification
- Automotive OEM Specifications

Special Features:

- Equipped with SBT2000 Windows based Automation Program.
- Excellent tool for analyzing low temperature pumpability of highly oxidized, 'used' and 'sooted' engine oils, lubricants and fuels.



05/4/8

Telephone: (989) 496-2309

FAX: (989) 496-3438

<http://www.savantgroup.com>

email: tannas@savantgroup.com

An ISO 9001 Certified Company